
Nestlé Health Science

Compat® Adapter

COMPAT® ADAPTER ENFit® SET → US FUNNEL TUBE

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EN **MD** This is a medical device. Adapter for connecting ENFit® administration sets to feeding tubes with legacy funnel ports. Intended for enterally fed patients. For use by healthcare professionals or by lay users and patients after proper training. ⚠ **Warnings:** For enteral use only. Ⓡ NOT for intravenous use. Do not use for more than 24 hours. Ⓢ Do not use if package is damaged. Do NOT use if labelling is incomplete or illegible. Always ensure appropriate patient monitoring. Small parts are a potential choking hazard if this device is not used as intended. Use under adult supervision. In case of drug administration to the patient, always verify that the medication delivered by the adapter is suitable for enteral feeding. Any reprocessing (especially sterilization) will compromise the structural integrity of the device. The device may become a source of infection and have altered performance after use and in case of reuse. This transition connector does not protect against inadvertent misconnection with connectors of non-enteral healthcare applications. **Important notes:** Ready to use. Use this transition connector ONLY if an enteral catheter equipped with an ENFit® connector is NOT available. DEHP: not intentionally added during the manufacturing process. Not formulated with BPA or natural rubber latex. Follow good clinical practice and all local and national regulations for transportation and disposal. Immediately report serious incident to the Manufacturer and to the Competent Authority in your country. Ⓢ This Instruction for Use is available at www.compat.com.

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C 00749-01 - 26MAY2021

Cedic srl

Via Liberazione 63/9, 20068 Peschiera Borromeo (MI), Italy
© Unless stated otherwise, all trademarks are used with permissions.
www.compat.com

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Nestlé Hong Kong Ltd.
Nestlé Singapore (Pte) Ltd.