
Compat® Adapter

COMPAT® ADAPTER ENFIT® TRANSITION KIT

تعليمات الاستخدام / Návod k použití / Brugsanvisning / Gebrauchsanweisung / Instructions for Use / Instrucciones de uso/ Käyttöohje / Mode d'emploi / Istruzioni per l'uso / Gebruiksaanwijzing / Bruksanvisning / Instrukcja użycia / Instruções de uso / Návod na použitie

REF 12334450
ANVISA 80062870007

Cedic srl
Via Liberazione 63/9,
IT-20068 Peschiera Borromeo (MI)

C00951-02
21SEP2023

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EN [MD] This is a medical device. **Adapter A: For connecting ENFit® administration sets to ENLock / Funnel feeding tubes. Adapter B: For connecting Luer / Oral / ENLock syringes to ENFit® medication ports.** Intended for enterally fed patients. For use by healthcare professionals or by lay users and patients after proper training. ⚠ **Warnings:** For enteral use only. Ⓢ Not for intravenous use. Do not use for more than 24 hours. Ⓢ Do not use if package is damaged. Do not use if labelling is incomplete or illegible. Always ensure appropriate patient monitoring. Small parts are a potential choking hazard if this device is not used as intended. Use under adult supervision. In case of drug administration to the patient, always verify that the medication delivered by the adapter is suitable for enteral feeding. Any reprocessing (especially sterilization) will compromise the structural integrity of the device. The device may become a source of infection and have altered performance after use and in case of reuse. This transition connector does not protect against inadvertent misconnection with connectors of non-enteral healthcare applications. **Important notes:** Ready to use. Use this transition connector only if a device equipped with an ENFit® connector is not available. DEHP: not intentionally added during the manufacturing process. Not formulated with BPA or natural rubber latex. Follow good clinical practice and all local and national regulations for transportation and disposal. Immediately report serious incident to the Manufacturer and to the Competent Authority in your country. 📄 This Instruction for Use is available at www.compat.com. To obtain further information please contact your sales representative or go to www.compat.com.

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