

# Compat ella<sup>®</sup>

## Compat Ella<sup>®</sup> Backpack

The Compat Ella<sup>®</sup> Backpack is intended for use with the Compat Ella<sup>®</sup> enteral feeding pump and the Compat Ella<sup>®</sup> Backpack Insert/Table Stand, by either adults or children.

### Features at a glance

#### Smart Design:


- Large zippered compartment makes it **easy to install and remove Compat Ella<sup>®</sup>**, administration set, and formula – without disconnection.
- Designed to work with the **backpack insert/table stand** for ultimate convenience.
- Administration set tubing exit on **left or right side**.
- **Padded back** for patient comfort with chest and hip straps for stability.

#### Convenient:

- Compatible with most tube feed containers **up to 1000mL**.
- Zippered flap for **quick access** to check or program Compat Ella<sup>®</sup>.
- Wide range of size adjustments makes it suitable for use by **children or adults**.
- Removable messenger strap gives **additional portability options** for wheelchair or as a messenger bag.

#### Robust and Modern:

- High quality materials, zippers and **water resistant** fabric protects Compat Ella<sup>®</sup> on the go.
- **Machine washable** at 30°C.
- Discreet **modern and compact design** in black and grey.




### Ordering Information

Product Name	Product Code	Qty/Box
Compat Ella <sup>®</sup> Backpack	12277736	1

Have you met Compat Ella<sup>®</sup>?

# Compat Ella® Backpack



## Components

- A Main compartment with dual zippers
- B 420D Nylon Ripstop fabric, with water resistant coating
- C Front flap for quick access to the pump
- D Zippered accessory compartment
- E Elastic webbing to secure small items
- F Removable wheelchair/messenger strap
- G External backpack hanger
- H 3D mesh padding for patient comfort
- I Adjustable shoulder straps
- J Adjustable chest and hips straps
- K Zippered side openings to allow tubing to exit from either left or right side
- L Side compartments for small accessories
- M Internal nutrition hanger

## Specifications

- Packing:** Individually packed in cases of 1 piece (40 x 24 x 4,5 cm).
- CE Mark:** Complies with the European Medical Device Directive 93/42 EEC.
- Legal manufacturer:** Axium MTech S.A., Allée de la Petite Prairie 4, CH-1260 Nyon, Switzerland.
- Ecological Footprint:** The 'Grüne Punkt' symbol signals that we contribute to the financing of recycling.
- Cleaning:** Machine wash at 30°C. Do not iron. Do not bleach. Do not tumble dry.
- Materials:** 420D Nylon Ripstop fabric.
- Country of origin:** Made in China.

To obtain further information please contact your Compat® representative or go to [www.compat.com](http://www.compat.com)  
®Unless stated otherwise, all trademarks are owned by Axium MTech S.A. or used with permission.

AXI00013 REV.0 14 May 2019

# Compat®