

# Compat ella®

## Compat Ella® Backpack Insert/Table Stand

The Compat Ella® Backpack Insert/Table Stand is a convenient dual-use pump and nutrition support intended for use either with the Compat Ella® Backpack for simple on-the-go use or as a stand-alone bed-side table stand.

### Features at a Glance

#### Convenient Dual-Use:


- Bed-side table stand: Small, portable, **convenient and low-profile**.
- Backpack insert: is easily transferred for use in the Compat Ella® Backpack without set disconnection, **for mobile use**.
- The extendable nutrition holder and adjustable straps makes it compatible with most tube feed containers up to 1500 mL.
- Easier to move around than a traditional IV pole.

#### Smart Design:

- **Easy access to all pump functions** when installed.
- Pump stays connectable to AC power
- Protects and holds Compat Ella® Pump securely.
- **"Kink-free" tubing management** with recessed tubing guides.
- 4 rubber foot pads **enhance stability** and reduce vibration.

#### Robust and Modern:

- Made of **durable polycarbonate**.
- **Discreet** white finish with rounded, curved design.




### Ordering Information

Product Name	Product Code	Qty/Box
Compat Ella® Backpack Insert/Table Stand	12277737	1

Have you met Compat Ella®?

# Compat Ella® Backpack Insert/Table Stand



## Components

- A Pump compartment
- B Nutrition compartment
- C Storage compartment for the AC adapter (when not connected to AC power)
- D Unique removable and adjustable stainless steel rod
- E Removable and adjustable nutrition container strap
- F "Kink-free" tubing management design with recessed tubing guides
- G Recessed design for AC adapter pump connection
- H Rubber foot pads

## Specifications

- Packing:** Individually packed in cases of 1 piece (17.5 × 30 × 17.5 cm).
- Storage:** Keep away from sunlight.
- CE Mark:** Complies with the European Medical Device Regulation 2017/745.
- Classification:** Class I.
- Legal manufacturer:** Axiom MTech S.A., Allée de la Petite Prairie 4, CH-1260 Nyon, Switzerland.
- Cleaning:** Hand wash.
- Materials:** Insert: Polycarbonate. Rod: Stainless steel. Strap: Polypropylene.
- Country of origin:** Made in China.

To obtain further information please contact your Compat® representative or go to [www.compat.com](http://www.compat.com)  
\*Unless stated otherwise, all trademarks are used with permission.  
This document does not replace the original Instructions for Use which must be read before use.

AXI00014 REV.3 05 DECEMBER 2022

# Compat®