

# Compat ella®

## Compat Ella® Backpack Insert/Table Stand

The Compat Ella® Backpack Insert/Table Stand is a convenient dual-use pump and nutrition support intended for use either with the Compat Ella® Backpack for simple on-the-go use or as a stand-alone bed-side table stand.

### Features at a glance

#### Convenient dual-use:


- **Bed-side table stand:** Small, portable, convenient and low-profile.
- **Backpack insert:** is easily transferred for use in the Compat Ella® backpack without set disconnection, for mobile use.
- The extendable nutrition holder and adjustable straps makes it compatible with most tube feed containers up to 1500mL.
- Easier to **move around** than a traditional IV pole.

#### Smart Design:

- **Easy access to all pump functions** when installed.
- Pump stays **connectable to AC power**
- **Protects and holds** Compat Ella® securely.
- “**Kink-free**” **tubing management** with recessed tubing guides.
- 4 rubber foot pads enhance **stability** and reduce vibration.

#### Robust and Modern:

- Made of **durable polycarbonate**.
- **Discreet** white finish with rounded, curved design.




### Ordering Information

Product Name	Product Code	Qty/Box
Compat Ella® Backpack/Table Stand	12277737	1

Have you met Compat Ella®?

# Compat Ella® Backpack Insert/Table Stand



## Components

- A Pump compartment
- B Nutrition Compartment
- C Storage compartment for the AC adapter (when not connected to AC power)
- D Unique removable and adjustable stainless steel rod
- E Removable and adjustable nutrition container strap
- F "kink-free" tubing management design with recessed tubing guides
- G Recessed design for AC adapter pump connection
- H Rubber foot pads

## Specifications

- Packing:** Individually packed in cases of 1 piece (17.5 x 30 x 17.5 cm).
- Storage:** Keep away from sunlight.
- CE Mark:** Complies with the European Medical Device Directive 93/42 EEC.
- Legal manufacturer:** Axium MTech S.A., Allée de la Petite Prairie 4, CH-1260 Nyon, Switzerland.
- Ecological Footprint:** The 'Grüne Punkt' symbol signals that we contribute to the financing of recycling.
- Cleaning:** Hand wash.
- Materials:** Insert: Polycarbonate. Rod: Stainless steel. Strap: Polypropylene.
- Country of origin:** Made in China.

To obtain further information please contact your Compat® representative or go to [www.compat.com](http://www.compat.com)  
\*Unless stated otherwise, all trademarks are owned by Axium MTech S.A. or used with permission.

AXI00014 REV.0 14 May 2019

# Compat®