

Compat ella®

Portable Enteral Feeding Pump

The Compat Ella® Pump is an innovative portable enteral feeding pump with a unique coloured user-friendly interface in local language.

Features at a Glance

- Continuous and interval feeding programs from 1 mL to 4000 mL, with or without dose limit
- Flow rates from 1 mL/h to 600 mL/h
- Rate and delivery progress is always visible in run mode

Easy to Learn and Understand

- A colour screen with easy to understand symbols for pump status, battery life, pump mode
- A clear U-shaped keypad with intuitive icon-based labelling
- On-screen full-text message line to convey pump status, help program nutritional therapies and **aid in problem solving**

Accurate and Safe

- Upstream and downstream occlusion alarms and air-in-line alarm
- +/- 5% delivery accuracy
- AAFF technology: Auto-Anti-Free-Flow to **promote patient safety**
- Separate delivery histories for 1) Current nutritional therapy and 2) Cumulative volume
- Keypad lock helps **avoid accidental changes**

Robust

- 5-year service interval
- Can be cleaned under warm running water

Mobile

- Lightweight and portable it can operate in any orientation
- 24-hour battery life at a delivery rate of 125 mL/h

Advanced Functions

- PDMS connection (Patient Data Management Systems)
- Interval mode delivers specific volumes at specified intervals
- Therapy lock to help **avoid accidental therapy changes**
- 30-day volume delivery history



Compat Ella® Pump Kit Contents

- Compat Ella® Pump
- AC power adapter
- Pole clamp
- Instructions for use

Compat Ella® Administration Sets

- Spike Set
- Combiset 1.5 L
- Universal Set
- Dual Spike Set

Compat Ella® Accessories

- Backpack
- Backpack Support Table Stand
- PDMS Connection Cable
- Nurse Call Connection Cable

Ordering Information

Product Name	Reference	Globe Code	Qty/Box
Compat Ella® Pump	12272020	12367032	1

Have you met Compat Ella®?

Compat Ella® Portable Enteral Feeding Pump

Pump Features

Delivery Accuracy: +/- 5%.

AAFF: Auto-Anti-Free-Flow (administration set based).

Keypad: Ergonomic U-shaped keypad with easy to identify icons.

Screen: Colour screen to highlight pump status.

Operating Languages: English, French, German, Italian, Spanish, Dutch, Swedish, Danish, Finnish, Norwegian, Czech, Slovakian, Portuguese and Polish.
User definable.

Continuous Mode: Continuous feed delivery. With or without total volume limit.

Flow Rate: 1-600 mL/h in 1 mL/h increments.

Volume Limit: 1-4000 mL in 1 mL increments.

Alarms: Air-in-line, separate upstream/downstream occlusion, door open, no cassette, end of dose, inactivity time-out, battery.

Volume This Therapy: Volume delivered for the currently active feeding nutritional therapy. User resettable.

Cumulative Volume: Total volume delivered since this parameter is reset by the user.

30-Day Volume History: Total volume delivered for each 24-hour period over the last 30 days.

Settings: Alarm volume, backlight intensity, backlight timeout, alarm at end of therapy (yes/no), keypad lock (on/off).

Priming: One touch automatic (17 mL) or manual.

Advanced Functions

Interval Mode: Delivery of specified bolus amounts in mL, at predefined intervals. Bolus from 1-2000 mL. Interval can be set from 15 minutes to 12 hours. With or without total volume limit.

Interval Mode Lockout: Interval mode is not available to the user.

Therapy Lock: Locks all the current nutritional therapy parameters so the user cannot change them.

Communication: PDMS or Nurse Call.

Date: User-adjustable date and time.

Patient Data Management System

Compat Ella® can connect to PDMS systems typically used in ICU - to allow tracking of feed delivery for individual patients in an automated manner. The physical connection is made via a custom serial cable. Specific software drivers are available upon request.

Warranty and Requirements

Warranty: 2 years.

Servicing: Only requires servicing every 5 years.

Servicing Agent: All service should be performed by authorized personnel only. Contact your Compat® representative or Compat® approved service provider.

Pump Specifications

Pump Weight: 499 g.

Battery: Li-Ion. Will run 24 hours at 125 mL/h with water. 6-Hour recharge time.

Pump Dimensions: L:140 mm × B:140 mm × H:62 mm.

Packaging: Individually packed in a box.

Storage: Store at room temperature. Keep away from sunlight and humidity.

Classification: Class IIa.

Mechanical and Fluid Ingress: IP25.

CE Mark: Complies with the European Medical Devices Directive 93/42/EEC. A CE Certification has been issued by an accredited Notified Body: TÜV SÜD.

Legal manufacturer: Axiom MTech S.A., Allée de la Petite Prairie 4, CH-1260 Nyon, Switzerland.

Country of origin: Made in Lithuania.

Compat Ella® Online

Website:

www.compat.com

Training Centre:

www.compatella.com

To obtain further information please contact your Compat® representative or go to www.compat.com

®Unless stated otherwise, all trademarks are used with permission.

This document does not replace the original Instructions for Use which must be read before use.

AXI00012 REV.3 05 DECEMBER 2022

Compat®