

# Compat ella®

## Portable Enteral Feeding Pump

The Compat Ella® Pump is an innovative portable enteral feeding pump with a unique coloured user-friendly interface in local language.

### Features at a glance

- Continuous and interval feeding programs from 1 mL to 4000 mL, with or without dose limit
- Flow rates from 1 mL/h to 600 mL/h
- Rate and delivery progress is always visible in run mode

#### Easy to learn and understand

- A colour screen with easy to understand symbols for pump status, battery life, pump mode
- A clear U-shaped keypad with intuitive icon-based labelling
- On-screen full-text message line to convey pump status, help program nutritional therapies and aid in problem solving

#### Accurate and safe

- Upstream and downstream occlusion alarms and air-in-line alarm
- +/-5% delivery accuracy
- AAFF technology: Auto-Anti-Free-Flow to enhance patient safety
- Separate delivery histories for 1) Current nutritional therapy and 2) Cumulative volume
- Keypad lock helps avoid accidental changes

#### Robust

- 5-year service interval\*
- Can be cleaned under warm running water

#### Mobile

- Lightweight and portable it can operate in any orientation
- 24-hour battery life at a delivery rate of 125 mL/h

#### Advanced Functions

- PDMS connection\*\* (Patient Data Management Systems)
- Interval mode delivers specific volumes at specified intervals
- Therapy lock to help avoid accidental therapy changes
- 30-day volume delivery history



#### Compat Ella® Pump Kit Contents

- Compat Ella® pump
- AC power adapter
- Pole clamp
- Instructions for use

#### Compat Ella® Administration Sets

- Spike set
- 1.5-Litre combiset
- Universal set
- Dual spike set

#### Compat Ella® Accessories

- Backpack
- Backpack support table stand
- PDMS connection cable\*\*
- Nurse call connection cable

### Ordering Information

Product Name	Reference	Globe Code	Qty/Box
Compat Ella® Pump	12272020	12367032	1

Have you met Compat Ella®?

# Compat Ella® Portable Enteral Feeding Pump

## Pump Features

**Delivery Accuracy:** +/-5%.

**AAFF:** Auto-Anti-Free-Flow (administration set based).

**Keypad:** Ergonomic U-shaped keypad with easy to identify icons.

**Screen:** Colour screen to highlight pump status.

**Operating Languages:** English, French, German, Italian, Spanish, Dutch, Swedish, Danish, Finnish, Norwegian, Czech, Slovakian, Portuguese and Polish.  
User definable.

**Continuous Mode:** Continuous feed delivery. With or without total volume limit.

**Flow Rate:** 1-600 mL/h in 1 mL increments.

**Volume Limit:** 1-4000 mL in 1 mL increments.

**Alarms:** Air-in-line, separate upstream/downstream occlusion, door open, no cassette, end of dose, inactivity time-out, battery.

**Volume This Therapy:** Volume delivered for the currently active feeding nutritional therapy. User resettable.

**Cumulative Volume:** Total volume delivered since this parameter is reset by the user.

**30-Day Volume History:** Total volume delivered for each 24-hour period over the last 30 days.

**Settings:** Alarm volume, backlight intensity, backlight timeout, alarm at end of therapy (yes/no), keypad lock (on/off).

**Priming:** One touch automatic (17 mL) or manual.

## Advanced Functions

**Interval Mode:** Delivery of specified bolus amounts in mL, at predefined intervals. Bolus from 1-2000 mL. Interval can be set from 15 minutes to 12 hours. With or without total volume limit.

**Interval Mode Lockout:** Interval mode is not available to the user.

**Therapy Lock:** Locks all the current nutritional therapy parameters so the user cannot change them.

**Communication:** PDMS\*\* or Nurse Call.

**Date:** User-adjustable date and time.

## \*\*Patient Data Management System

Compat Ella® can connect to PDMS systems typically used in ICU - to allow tracking of feed delivery for individual patients in an automated manner. The physical connection is made via a custom serial cable. Specific software drivers are available upon request.\*

## Warranty and Requirements

**Warranty:** 2 years.

**Servicing:** Only requires servicing every 5 years.\*

**Servicing Agent:** All service should be performed by authorized personnel only. Contact your Compat® representative or Compat® approved service provider.

## Pump Specifications

**Pump Weight:** 499 g.

**Battery:** Li-Ion. Will run 24 hours at 125 mL/h with water. 6-Hour recharge time.

**Pump Dimensions:** L:140 mm x B:140 mm x H:62 mm.

**Pump Kit Packaging:** Individually packed in a box.

**Storage:** Store at room temperature. Keep away from sunlight and humidity.

**Device Classification:** Class IIa.

**Mechanical and Fluid Ingress:** IP25.

**CE Mark:** Complies with the European Medical Device Directive 93/42 EEC. A CE Certification has been issued by an accredited Notified Body: TÜV SÜD.

**Legal manufacturer:** Axiom MTech S.A., Allée de la Petite Prairie 4, CH-1260 Nyon, Switzerland.

**Ecological Footprint:** The 'Grüne Punkt' symbol signals that we contribute to the financing of recycling.

**Country of origin:** Made in Lithuania.

## Compat Ella On-line

**Website:**

[www.compat.com](http://www.compat.com)

**Training Centre:**

[www.trainingcentre.compatella.com](http://www.trainingcentre.compatella.com)

To obtain further information please contact your Compat® representative or go to [www.compat.com](http://www.compat.com)  
\*Unless stated otherwise, all trademarks are owned by Axiom MTech S.A. or used with permission.

AXI00012 REV.0 14 May 2019

# Compat®