

# Compat<sup>®</sup> Adapter

## Compat<sup>®</sup> Adapter ENFit<sup>®</sup> Set > US Funnel Tube

The Compat<sup>®</sup> Adapter ENFit<sup>®</sup> Set > US Funnel Tube is intended to connect an ENFit<sup>®</sup> administration set with a funnel tube.

### Features at a Glance

The Compat<sup>®</sup> Adapter ENFit<sup>®</sup> Set > US Funnel Tube is placed between the administration set and the feeding tube.

**Ease of Use:** The Compat<sup>®</sup> Adapter ENFit<sup>®</sup> Set > US Funnel Tube is designed to **facilitate the transition to ENFit<sup>®</sup>** by allowing compatibility between ENFit<sup>®</sup> administration sets and feeding tubes with legacy funnel ports.

**Safer Connectors:** The Compat<sup>®</sup> Adapter ENFit<sup>®</sup> Set > US Funnel Tube includes an ENFit<sup>®</sup> connector made according to ISO 80369-3. ENFit<sup>®</sup> is designed to be IV incompatible in order to help **promote patient safety** by helping to avoid accidental misconnections.

**Safe Materials:** The Compat<sup>®</sup> Adapter ENFit<sup>®</sup> Set > US Funnel Tube is **DEHP, BPA and Latex-Free\***.

**Compatibility:**

- Administration Sets: ENFit<sup>®</sup>
- Enteral feeding tubes: Funnel\*\*




### Ordering Information

Product Name	Product Code	Qty/Box
Compat <sup>®</sup> Adapter ENFit <sup>®</sup> Set > US Funnel Tube	12298016	30 pcs

*Patient Care. Our Priority.*

# Compat® Adapter ENFit® Set > US Funnel Tube



## Components

**A** Adapter ENFit® Set > US Funnel Tube

## Specifications

<b>Sterilization:</b>	Non-sterile product.
<b>Packing:</b>	Packed in blisters of 5 pieces. 30 pieces per case (11,5 × 9,5 × 9 cm).
<b>Storage:</b>	Store at room temperature. Keep away from sunlight and humidity.
<b>Shelf life:</b>	3 years from production date. Expiry date on the package.
<b>Duration of use:</b>	24 hours. Single Use.
<b>ISO Standard:</b>	ENFit® connector made according to ISO 80369-3.
<b>CE Mark:</b>	Complies with the European Medical Device Regulation 2017/745.
<b>Classification:</b>	Class I.
<b>Legal manufacturer:</b>	Cedic Srl, Via Liberazione 63/9, IT-20068 Peschiera Borromeo (MI).
<b>Sustainability:</b>	Environmental management system complies with ISO 14001.
<b>Materials:</b>	All materials DEHP, BPA and Latex-free*. Additional details on request.
<b>Country of origin:</b>	Made in Italy.

\*DEHP: Not intentionally added during the manufacturing process. Not formulated with BPA or natural rubber latex.

\*\*Some incompatibilities may exist.

To obtain further information please contact your Compat® representative or go to [www.compat.com](http://www.compat.com)

\*Unless stated otherwise, all trademarks are used with permission.

This document does not replace the original Instructions for Use which must be read before use.

AXI01702 Rev1 14 MARCH 2024

# Compat®