

Compat[®] Adapter

Compat[®] Adapter ENFit[®] to Oral/Luer Syringe

The Compat[®] Adapter ENFit[®] to Oral/Luer Syringe is intended for connecting an enteral Oral, Luer or ENLock syringe to an enteral giving set equipped with an ENFit[®] medication port. It is also intended for use with a feeding tube with ENFit[®] connector.

Features at a Glance

The Compat[®] Adapter ENFit[®] to Oral/Luer Syringe is placed between the feeding tube and the syringe or between the medication port and the syringe.


Ease of Use: The Compat[®] Adapter ENFit[®] to Oral/Luer Syringe is designed to **facilitate the transition to ENFit[®]** by allowing compatibility between ENFit[®] medication ports and feeding tubes with legacy syringes.

Safer Connectors: The Compat[®] Adapter ENFit[®] to Oral/Luer Syringe features an ENFit[®] connection in accordance with ISO 80369-3. ENFit[®] is designed to be IV incompatible in order to help **promote patient safety** by helping to avoid accidental misconnections.

Safe Materials: The Compat[®] Adapter ENFit[®] to Oral/Luer Syringe is **sterile, DEHP, BPA and Latex-Free***.

Compatibility:

- Enteral feeding tubes: ENFit[®]
- Medication port: ENFit[®]
- Syringes: ENLock, Oral, Luer




Ordering Information

Product Name	Product Code	Qty/Box
Compat [®] Adapter ENFit [®] to Oral/Luer Syringe	12459535	30 pcs

Patient Care. Our Priority.

Compat® Adapter ENFit® to Oral/Luer Syringe



Components

- A Oral/Luer adapter
- B Removable cap

Specifications

Sterilization:	Sterilized with Ethylene Oxide (EO).
Packing:	Individually packed in cases of 30 pcs (11,5 x 9,5 x 9 cm).
Storage:	Store at room temperature. Keep away from sunlight and humidity.
Shelf life:	3 years from production date. Expiry date on the package.
Duration of use:	24 hours. Single Use.
ISO Standard:	ENFit® connector made according to ISO 80369-3.
CE Mark:	Complies with the European Medical Device Directive 93/42 EEC. A CE Certification has been issued by an accredited Notified Body: TÜV SÜD (CE0123).
Classification:	Class Is.
Legal manufacturer:	Cedic Srl, Via Liberazione 63/9, IT-20068 Peschiera Borromeo (MI).
Sustainability:	Environmental management system complies with ISO 14001:2015.
Materials:	All materials DEHP, BPA and Latex-free*. Additional details on request.
Country of origin:	Made in Italy.

*DEHP: Not intentionally added during the manufacturing process. Not formulated with BPA or natural rubber latex.

To obtain further information please contact your Compat® representative or go to www.compat.com

®Unless stated otherwise, all trademarks are used with permission.

This document does not replace the original Instructions for Use which must be read before use.

AXI01535 REV0 22 NOVEMBER 2022

Compat®