

Compat[®] ENject

Enteral Syringes with ENFit[®] Connector

The Compat[®] ENject Enteral Syringes with ENFit[®] connector are indicated for use with ENFit[®] enteral access medical devices.

Features at a Glance

The Compat[®] ENject Syringes feature a transparent body with graduated markings, a purple plunger, an ENFit[®] connector and a cap.

Convenience: Compat[®] ENject Syringes are available in 8 sizes from 1 mL to 100 mL. Low Dose Tip (LDT) syringes are available in sizes 1 mL, 2.5-3 mL and 5 mL. The plunger **slides easily in the syringe body** and the syringe features a **slim body for ease of handling**. The **cap** allows users to **prepare feeds or medication in advance** and **prevents syringe from leaking**.


Patient Safety: Compat[®] ENject Syringes feature easy-to-read volume markings to **facilitate accurate delivery of fluids**. LDT syringes feature a low dead space reverse ENFit[®] connector to **help ensure accurate dosing** when delivering small medication volumes. All Compat[®] ENject Syringes feature ENFit[®] connectors designed to be IV incompatible in order to **help promote patient safety** and avoid accidental misconnections.

Safe Materials: Compat[®] ENject syringes are sterile, DEHP, BPA and Latex-Free¹.

Compatibility:

Enteral feeding tubes: ENFit[®]

Medication administration ports: ENFit[®]




Ordering Information

Product Name	Product Reference	Units per Box
Compat [®] ENject Enteral Syringe 1 mL, Low Dose Tip	P00103	100
Compat [®] ENject Enteral Syringe 2.5-3 mL, Low Dose Tip	P00104	100
Compat [®] ENject Enteral Syringe 5 mL, Low Dose Tip	P00105	100
Compat [®] ENject Enteral Syringe 10 mL	P00106	100
Compat [®] ENject Enteral Syringe 20 mL	P00107	50
Compat [®] ENject Enteral Syringe 30 mL	P00198	50
Compat [®] ENject Enteral Syringe 50-60 mL	P00108	30
Compat [®] ENject Enteral Syringe 100 mL	P00199	30

Patient Care. Our Priority.

Enteral Syringes with ENFit® Connector



Components

- A Transparent body with graduations
- B Purple plunger
- C 1 mL – 30 mL, 100 mL syringes: centric ENFit® connector
- D 50-60 mL syringe: eccentric ENFit® connector
- E ENFit® cap

Specifications

- Sterility:** Sterilized with Ethylene Oxide (EO).
- Packing:** Individually packed in cases of 30 to 100 pieces.
- Storage:** Store dry at room temperature and away from sunlight.
- Shelf life:** 5 years from production date. Expiry date on the package.
- Duration of use:** For single use.
- Size available:** From 1 mL to 100 mL. See table on page 1.
- ISO Standard:** ENFit® connector made according to ISO 80369-3.
- CE Mark:** Complies with the European Medical Devices Directive 93/42/EEC. A CE certification has been issued by an accredited Notified Body: TÜV Rheinland (CE0197).
- Classification:** Class Is.
- Legal manufacturer:** Jiangsu Micsafe Medical Technology Co., Ltd Xituan Industrial Park, Dafeng District, Yancheng, 224125, Jiangsu, China.
- EC representative:** Share Info Consultant Service LLC Repräsentanzbüro, Heerdter Lohweg 83, 40549 Düsseldorf.
- Materials:** Polypropylene. All materials DEHP, BPA and Latex-free¹. Additional details on request.
- Country of origin:** Made in China.

¹DEHP and BPA are not intentionally added during the manufacturing process. Not formulated with natural rubber latex.

To obtain further information please contact your Compat® representative or go to www.compat.com

* Unless stated otherwise, all trademarks are used with permission.