

Compat[®] Adapter

Compat[®] Adapter ENFit[®] Transition Kit

The Compat[®] Adapter ENFit[®] Transition Kit is intended for allowing the connection of ENFit[®] administration sets to legacy feeding tubes and syringes.

Features at a Glance

The Compat[®] Adapter ENFit[®] Transition Kit features an ENLock/Stepped Adapter (adapter A) and a Luer/Oral Adapter (adapter B). It is designed to **allow compatibility between ENFit[®] administration sets and most enteral feeding tubes and syringes.**

Food Safe Materials: The Compat[®] Adapter ENFit[®] Transition Kit is DEHP, BPA and Latex-Free*.

Compatibility:


Adapter A: ENLock/Stepped Adapter

- Administration sets: ENFit[®]
- Enteral feeding tubes: ENLock, Funnel**

Adapter B: Luer/Oral Adapter

- Syringes: Luer, Oral, ENLock
- Enteral feeding tubes: ENFit[®]
- Medication port: ENFit[®]

Adapter A:



Adapter B:




Ordering Information

Product Name	Product Code	Qty/Box
Compat [®] Adapter ENFit [®] Transition Kit	12334450	10 pcs


Patient Care. Our Priority.

Compat® Adapter ENFit® Transition Kit

Adapter A:



Adapter B:



Components

Adapter A:

- I ENLock/Stepped Adapter
- II Removable cap

Adapter B:

- I Luer/Oral Adapter
- II Removable cap

Specifications

Sterile:	Non-sterile product.
Packing:	Individually packed in cases of 10 pieces (11,5 x 9,5 x 9 cm).
Storage:	Store at room temperature. Keep away from sunlight and humidity.
Shelf life:	3 years from production date. Expiry date on the package.
Duration of use:	24 hours. Single Use.
CE Mark:	Complies with the European Medical Device Regulation 2017/745.
Classification:	Class I.
Legal manufacturer:	Cedic Srl, Via Liberazione 63/9, IT-20068 Peschiera Borromeo (MI).
Sustainability:	Environmental management system complies with ISO 14001:2015.
Materials:	All materials DEHP, BPA and Latex-free*. Additional details on request.
Country of origin:	Made in Italy.

*DEHP and BPA are not intentionally added during the manufacturing process. Not formulated with natural rubber latex.

**Some incompatibilities may exist.

This document does not replace the original Instructions for Use which must be read before use. To obtain further information please contact your Compat® representative or go to www.compat.com

*Unless stated otherwise, all trademarks are used with permission.

AXI01533 REV.0 07 NOVEMBER 2022

Compat®